

Brain Tumor Classification Based on Deep CNN using MRI images

Capstone project

Momen Albawarid (21110399)

Supervisor: Eng. Ayah Karajeh

From the CRP to The Capstone Project

Why the Medical Field?

The Problem Statement

1.

The Need for Automated Clinical Support

2.

The "Black Box" Problem in Medical AI

3.

Model Generalization and Stability

4.

Technical and Resource Constraints

Research Phase (CRP)

What We Learned

1.

Architecture Efficiency vs. Performance

The study evaluated five pre-trained models

Xception ResNet50 Inception V3 VGG16 **MobileNet**

on an extensive MRI dataset

MobileNet emerged as the most efficient model

**balancing high accuracy with
the lowest memory footprint**

2.

The Negative Impact of Data Augmentation

3.

**Image Enhancement
is Key**

Capstone Phase

What Did I do

High Accuracy

Optimizing MobileNetV2 to achieve >90% accuracy

Explainable AI

Integrate Score-CAM for clean, localized tumor
heatmaps




Efficiency

Ensure CPU inference under 1,5s for rapid triage.

Usability

Desktop GUI for users

System Architecture



Application Walkthrough

Brain Tumor Classification

Select an MRI image, analyze, then review Score-CAM localization.

Input MRI

Click 'Browse MRI' to select an image

No file selected

Prediction: —

Confidence: —

Visualization

Original


Score-CAM Overlay

Browse MRI

Analyze

New Image

Ready



Brain Tumor Classification

Select an MRI image, analyze, then review Score-CAM localization.

Input MRI

Click 'Browse MRI' to select an image


C:/Users/DELL/PycharmProjects/CapstoneProject2/test_image_meningioma.jpg

Prediction: —

Confidence: —

Visualization

Original



Score-CAM Overlay

—

Browse MRI

Analyze

New Image

Image loaded. Press Analyze.

Brain Tumor Classification

Select an MRI image, analyze, then review Score-CAM localization.

Input MRI

Click 'Browse MRI' to select an image

C:/Users/DELL/PycharmProjects/CapstoneProject2/test_image_meningioma.jpg

Prediction: —

Confidence: —

Visualization

Original

Score-CAM Overlay

Please wait

Analyzing MRI...

Running Score-CAM on CPU...

Browse MRI

Analyze

New Image

Analyzing...

Brain Tumor Classification

Select an MRI image, analyze, then review Score-CAM localization.

Input MRI

Click 'Browse MRI' to select an image


C:/Users/DELL/PycharmProjects/CapstoneProject2/test_image_meningioma.jpg

Prediction: Meningioma Tumor


Confidence: 1.000

Visualization

Original



Score-CAM Overlay



Browse MRI

Analyze

New Image

Done. Click New Image to analyze another.

Brain Tumor Classification

Select an MRI image, analyze, then review Score-CAM localization.

Input MRI

Click 'Browse MRI' to select an image

No file selected

Prediction: —

Confidence: —

Visualization

Original

Score-CAM Overlay

Browse MRI

Analyze

New Image

Ready

Performance and Results Evaluation

The final model achieved a 93% Test Accuracy
on a hold-out set of 3,961 images

	Precision	Recall	F1-Score
Glioma	0.94	0.96	0.95
Meningioma	0.9	0.86	0.88
No Tumor	0.92	0.95	0.93
Pituitary	0.96	0.98	0.97
Weighted Avg	0.93	0.93	0.93

Resource & Speed Efficiency

A key goal was making the tool run on standard consumer hardware

Inference Speed

The average analysis takes less than 1.5 seconds on a standard CPU.

Memory Footprint

The model uses only ~200MB of RAM, compared to ~1400MB required by heavier models like Xception.

Startup Optimization

The system is designed to load the .h5 model once at startup, which removes any "lag" during subsequent image checks.

Software Reliability (Manual Testing)

A key goal was making the tool run on standard consumer hardware

Test ID	Scenario	Action Performed	Expected Result	Status
1	Standard Workflow	Load a valid .jpg MRI image and click "Analyze".	Prediction and Heatmap appear side-by-side. Status updates to "Done".	PASS
2	File Compatibility	Load a valid .png MRI image.	Image loads and analyzes correctly (same as JPG).	PASS
3	Invalid File Type	Attempt to load a .txt or .pdf file.	File dialog automatically filters view to only show supported formats (.jpg, .png), preventing invalid selection.	PASS
4	Concurrency	Click "Analyze", then immediately drag the window.	Window moves smoothly; no "Not Responding" freeze.	PASS
5	Empty State	Click "Analyze" without selecting an image first.	The 'Analyze' button remains disabled (grayed out) until a valid image is loaded.	PASS
6	Rapid Clicks	Click "Analyze" button 5 times quickly.	Button disables after 1st click; only one analysis runs.	PASS
7	Session Reset	Analyze Image A, then immediately load Image B.	Previous results (prediction/heatmap) clear; new image loads ready for analysis.	PASS

Challenges and Limitation

The "Grid" Issue in Explainable AI (XAI)

Resource and Computing Constraints

Data Quality and Diversity

Future Work

Clinical & Machine Integration

Direct MRI Integration

3D Volumetric Analysis

Automated Medical Reporting

Technical Advancements

Segmentation Techniques

Advanced Architectures

Thank you for listening